




PathAssist™
LIGHT SEEKER™

The Entellus™ PathAssist™
Light Seeker™ is a
versatile, reusable,
confirmation tool
to locate the
frontal sinus.



Lighted confirmation
of frontal sinus

An easy-to-use tool that combines the tactile feedback of
a sinus seeker with the confidence of illuminated guidance.

- Confirms correct entry of the frontal recess and location of the frontal sinus
- Provides confirmation during in-office procedures
- Simple, economical alternative to costly guidance systems for OR FESS procedures
- Reusable with convenient multiple light source connections

Angle: 78 degrees
Ball Tip: 1.5mm

A practical, easy-to-use, confirmation tool for sinus procedures in the office or OR.

A new addition to Entellus™ PathAssist™ tools, the Entellus Light Seeker is a cost-effective, reusable instrument that enables ENTs to confidently navigate nasal pathways – without the need for costly additional guidance technologies.

The Light Seeker illuminates the correct entry to the frontal recess while providing transcutaneous confirmation of the frontal sinus.

An ideal general-use tool for in-office or OR sinus procedures.

Whether used to identify the frontal sinus prior to irrigation or for endoscopic exploration, the Light Seeker is a versatile, economical tool that can supplement or replace costly visualization technologies.

The Light Seeker has a standard ACMI connection. Two converter adapters are provided for added connectivity to existing light cables.

To learn more about the Light Seeker visit www.entellusmedical.com

About Entellus Medical

Entellus Medical delivers innovative, high quality minimally invasive therapeutic solutions to healthcare providers and their patients who suffer from sinusitis.



Light Seeker Indications for Use:

The PathAssist Light Seeker is intended to locate, illuminate within, and transilluminate across nasal and sinus structures, including the frontal, ethmoid and maxillary sinuses, in adults aged 18 and over. Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and adverse events. Refer to PathAssist Light Seeker Instructions for Use (IFU) for complete instructions.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.