

Compare XprESS™ Multi-Sinus Dilation System

PRODUCT (Download Info)	XprESS Ultra Super sleek profile delivers optimal visualization and enables navigation of tight anatomy during both office and OR cases Learn More	XprESS LoProfile Broader range of balloon sizes allows physician to address special case patient needs Learn More	XprESS Pro Wider lumen for enhanced integrated suction/irrigation Learn More
Ball Tip	1.5mm	1.75mm	2.0mm
Diameter of hypotube at 1cm	1.0mm	1.2mm	1.5mm
Irrigation	yes	yes	yes
Suction comparable to	2.5 French	4.0 French	5.0 French
Compatible with CT Image Guidance	NA	Yes	Yes
Available standard balloon sizes	5 x 20mm 6 x 20mm	5 x 20mm 6 x 20mm 7 x 20mm	6 x 18mm 7 x 18mm
Additional balloon sizes	5 x 8mm 6 x 8mm	5 x 8mm 6 x 8mm	6 x 8mm
Integrated PathAssist LED Light Fiber	Included	Included	Add-on
Bend Options	Frontal Sphenoid Maxillary 135 Maxillary 120 Eustachian tube	Frontal Sphenoid Maxillary 135 Maxillary 120 Eustachian tube	Frontal Sphenoid Maxillary 135 Maxillary 120 Eustachian tube
Intuitive inflation syringe	Included	Included	Included
Single Box Packaging	yes	yes	yes

Order Details

To order, please call 866-620-7615 or [contact us](#).

CATALOG NO.	DESCRIPTION
ULF-105	XprESS™ Ultra 5 mm x 20 mm
ULF-106	XprESS™ Ultra 6 mm x 20 mm
ULF-205	XprESS™ Ultra 5 mm x 8 mm
ULF-206	XprESS™ Ultra 6 mm x 8 mm
LPLF-105-I	XprESS™ LoProfile 5 mm x 20 mm
LPLF-106-I	XprESS™ LoProfile 6 mm x 20 mm
LPLF-107-I	XprESS™ LoProfile 7 mm x 20 mm
LPLF-205-I	XprESS™ LoProfile 5 mm x 8 mm
LPLF-206-I	XprESS™ LoProfile 6 mm x 8 mm
JD-106	XprESS™ Pro 6 mm x 18 mm
JD-107	XprESS™ Pro 7 mm x 18 mm
JD-206	XprESS™ Pro 6 mm x 8 mm