

XprESS™ ENT Dilation System Technical Specifications

PRODUCT (Download)	XprESS LoProfile Broad range of balloon sizes allows physician to address special case patient needs
Ball Tip	1.75mm
Diameter of hypotube at 1cm	1.2mm
Irrigation	yes
Suction comparable to	4.0 French
Compatible with CT Image Guidance	NA
Available standard balloon sizes	5 x 20mm 6 x 20mm 7 x 20mm
Additional balloon sizes	5 x 8mm 6 x 8mm
Integrated LED Light Fiber	Included
Bend Options	Frontal Sphenoid Maxillary 135 Maxillary 120 Eustachian tube
Intuitive inflation syringe	Included
Single Box Packaging	yes

Order Details

To order, please call 866-620-7615 or [contact us](#).

CATALOG NO.	DESCRIPTION
LPLF-105-I	XprESS™ LoProfile 5 mm x 20 mm
LPLF-106-I	XprESS™ LoProfile 6 mm x 20 mm
LPLF-107-I	XprESS™ LoProfile 7 mm x 20 mm
LPLF-205-I	XprESS™ LoProfile 5 mm x 8 mm
LPLF-206-I	XprESS™ LoProfile 6 mm x 8 mm